

Fresenius 2008T Bibag Door Error – Quick Reference Guide

■ Symptom

Display Message: 'bibag: Door Error'

Occurs during initialization or operation when the door sensor signal is not detected as closed.

■ Root Causes

1. Door not fully latched – gray latch not pushing the microswitch completely.
2. Broken/misaligned latch tab – plastic tab on the door worn or cracked.
3. Faulty microswitch – electrical failure or mechanical sticking.
4. Loose/corroded harness connection – wiring between door switch and bibag board.
5. Bibag interface board input fault – rare logic fault on control input.

■ Labeled Diagram Reference

[A] Gray Latch Tab – locks the door and presses the switch.

[B] Microswitch Button – round gray/white button beneath latch; sends signal.

[C] O-Ring Ports – bicarbonate bag seals (top left/right).

[D] Bibag Door Hinge – where misalignment may occur.

[E] Harness Connector – wiring to bibag board.

■ Inspection Steps

1. Visual Check: Ensure latch [A] fully presses switch [B].
2. Manual Actuation: Press microswitch button [B] – should click.
3. Continuity Test: Door open → OL; Door closed → 0Ω.
4. Inspect Harness [E]: Check for corrosion or loose pins.
5. Jumper Test: Temporarily short connector to simulate closed door.
6. Input Verification: Biomed Mode → Diagnostics → Input Test → Verify toggle between 0 & 1.

■ Corrective Actions

| Issue | Corrective Action |

| Door latch cracked/worn | Replace Bibag Door Assembly |

| Microswitch defective | Replace Door Switch |

| Corroded harness | Clean/replace harness |

| Bibag board fault | Replace Bibag Interface Board |

■ Technician Notes

- Microswitch is behind latch [A]; button [B] visible under latch.
- No click = mechanical failure.
- Avoid forcing door; can damage hinge or tab.
- Dry moisture before testing.